

How to use this table:

- 1) Look up the row with the chemical you are using
- 2) follow the row to find a glove without red background
- 3) Check the BT time
- 4) If that time is higher than the intended time you are going to use it, that glove is suitable
- 5) if not, continue to the next non-red glove

Example - Formaldehyde 37%:

The best is Microflex (Level 6), if you are going to handle formaldehyde less than 87 min you can use the orange or green. Normally you should replace gloves as soon as you know they are contaminated with something that might penetrate it, so a not very good glove may be OK.

If you are cleaning up spill you need a very good glove

BTT = Breakthrough Time (min)

Chem protection			Glove	SHIELDSkin™ Orange Nitrile™ 260		Microflex 93-850		BTT = Breakthrough Time (min)	
	time	level		EN ISO 374-5:2016 Level 3 VIRUS	EN ISO 374-1:2016 TYPE B JKPT	EN ISO 374-5:2016 VIRUS	EN ISO 374-1:2016 TYPE B JKOPT	Level Abena Excellent	BTT Abena Excellent
Not good	0	0	Virus protection	Level Orange Nitrile	BTT Orange Nitrile	Level Microflex 850	BTT Microflex 50	Level Abena Excellent	BTT Abena Excellent
Short time use only	10	1							
	30	2	Chemical protection	Level Orange Nitrile	BTT Orange Nitrile	Level Microflex 850	BTT Microflex 50	Level Abena Excellent	BTT Abena Excellent
Acceptable	60	3							
	120	4	Chemical protection	Level Orange Nitrile	BTT Orange Nitrile	Level Microflex 850	BTT Microflex 50	Level Abena Excellent	BTT Abena Excellent
Good	240	5							
Very good	480	6	Chemical protection	Level Orange Nitrile	BTT Orange Nitrile	Level Microflex 850	BTT Microflex 50	Level Abena Excellent	BTT Abena Excellent
For a list of additional tested chemicals, see other list or product data sheet									
Chemical	EN 374 1 Code	CAS	Level Orange Nitrile	BTT Orange Nitrile	Level Microflex 850	BTT Microflex 50	Level Abena Excellent	BTT Abena Excellent	
Acetic Acid 10%	64-19-7		0	0	0	7	0	8	Not tested
Acetonitrile	C 75-05-8		Not tested		0	1	0	1	Not tested
Ammonium Hydroxide 25%	O 1336-21-6		0	0	0	0	0	2	0
Dichloromethane	D 75-09-2		6	480	6	480	Not tested		6
Diethyl ether	60-29-7		1	16	2	31	3	101	1
Dimethyl Sulfoxide 99% (DMSO)	67-68-5		0	0	0	0	0	1	0
Ethanol 99%	64-17-5		0	0	0	0	0	1	0
Formaldehyde 10%	<T 50-00-0		0	0	0	1	0	3	0
Formic acid	64-18-6		0	0	0	4	Not tested		Not tested
n-Hexane	110-54-3		2		3	91	6	480	Not tested
Hydrogen Peroxide 12%	<P 7722-84-1		1	28	3	97	6	480	6
Isopropanol	67-63-0		4		4	130	4	204	Not tested
2-Mercaptoethanol	60-24-2		Not tested		6	480	Not tested		Not tested
Methanol 20%	<A 67-56-1		6	480	6	480	6	480	Not tested
Phenol 4%	108-95-2		2		2	46	5	242	Not tested
Rely+On Virkon	70693-62-8		3		3	70	Not tested		Not tested
Sodium Hypochlorite 13%	X 7681-52-9		0	1	0	1	Not tested		Not tested
Sulphuric Acid 50%	<L 7664-93-9		0	0	0	5	0	6	Not tested
Tetrahydrofuran	H 109-99-9		3		Not tested		Not tested		Not tested
Triethylamine	121-44-8		0	4	0	5	0	5	Not tested
o-Xylene	95-47-6		Not tested		1	22	Not tested		Not tested
			6	480	6	480	Not tested		Not tested
			Not tested		Not tested		Not tested		6
			6	480	6	480	6	480	6
			6	480	6	480	Not tested		Not tested
			6	480	6	480	Not tested		Not tested
			6	480	6	480	6	480	Not tested
			0	0	0	8	1	12	0
			0	0	0	0	0	2	Not tested
			Not tested		0	1	0	3	0
			Not tested		1	23	Not tested		Not tested
			0	0	0	1	Not tested		Not tested
			Not tested		Not tested		Not tested		0
			0		0	4	0	5	Not tested